


# The Advantages of Breasy



Breasy's revolutionary design combines the inhalation mask with Breasy's supplied pacifier or the child's own pacifier, thus facilitating familiarity and comfort. It can be used in various scenarios such as in the E.R. for sedation or ventilation purposes, MRI procedures and of course for inhalation procedures and the treatment of various breathing illnesses.



Concurrently Breasy targets home use, where children are encouraged to play with the mask prior to treatment and make it a part of their everyday life. When the need arises there is no resistance to applying the mask and combined with the calming effect of the pacifier the child is calm and resistance to treatment is dramatically reduced.



Breasy is suitable for use with all standard nebulizers and pacifiers currently in market.

## The Science behind Breasy's design

The suction effect caused by using the pacifier creates a perfect fit seal around the child's mouth and nose, allowing for efficient airflow from the holding chamber to the child's lungs. This perfect seal fit is also what allows Breasy to be the only mask that requires no holding straps during treatment.

Breasy unfolds an innovative combination: a significantly lower resistance to treatment along with a perfect fit design, leading to substantially improved treatment efficiency when it matters the most.

